

AUTO REF/KERATOMETER
GR-3500KA
GR-3300K

| Specification | |
|-------------------------------|---|
| Model | GR-3500KA GR-3300K |
| Refraction Measurement | Sphere -22~+30D (VD=0mm) |
| | -30~+22D (VD=12mm) (0.01/0.12/0.25DStep) |
| | Cylinder 0~±10D (0.01/0.12/0.25DStep) |
| | Axis angle 1~180° (1°Step) |
| Measurement of Corneal Radius | Corneal radius 5.0~10.0mm (0.01mmStep) |
| | Refractive power 33.75~67.5D (Corneal Refractive Index n=1.3375) |
| | (0.01/0.12/0.25DStep) |
| | Cylindrical power 0~±10D Axis angle 1~180° (1°Step) |
| Corneal Measured Area | φ2.8mm (Ring Measurement/at 8.0mm of Corneal Radius) |
| | φ7.0mm (4 Point Peripheral Measurement/at 8.0mm of Corneal Radius) |
| Pupil Diameter Measurement | φ2~φ8mm (0.1mm step) |
| Vertex Distance | 0, 10, 12, 13.5, 15mm |
| Minimum Pupil Diameter | φ2.2mm |
| Pupillary Distance | Measurement Range 0~85mm (1mmStep) |
| Auto Start | ○ |
| Auto Focus/Auto Track | ○ × |
| Printer | Thermal printer with Automatic Cutter (Width 57mm) |
| Internal Monitor | 5.7 inch LCD Display (Color) |
| Movable Distance | Back/Force ±17mm Right/Left ±43mm Up/Down ±17mm |
| Movable Distance of Chinrest | ±30mm |
| Overall Dimension | (W) 260mm × (D) 465mm × (H) 453mm |
| Weight | About 20kg About 17kg |
| Output | RS-232C, USB2.0 Interface |
| Rated Voltage | 100~240V 50/60Hz |
| Consumption | 90VA 80VA |
| Power Save | OFF, 3, 5, 10 min. (Selectable) |


Accessories

- Test eye / 1 pc ● Power cable / 1 pc ● Printer paper / 3 rolls
- Fuse 2A / 2 pcs ● Chinrest paper / 1 pack (1000 sheets)
- Chinrest pin / 2 pcs ● Dust cover / 1 pc ● Contact lens holder / 1 pc
- Operation Manual / 1 pc ● USB Driver CD / 1 pc

Option

- MDC-1 The Measurement Data Collection software
- Electric table ● Chair

Overall Dimension



- Specifications and design are subject to change without prior notice for improvement.
- The screen is composed photograph.
- The color of instrument on the catalogue and the real product might be different.
- We judge that the LCD Monitor is qualified if the total "lit pixel" and/or "missing pixel" is less than five (excluding non-sharp one and less than half one).
- Production and Distribution Registration Number GR-3500KA: 37B2X10001000037
GR-3300K : 37B2X10001000039

- Specifications and design are subject to change without prior notice for improvement.
 - The color of instrument on the catalogue and the real product might be different.
- CAUTION** For the proper and safety use, please read the Instruction Manual thoroughly before using the product.

● Agency

Rexxam CE 0197
Rexxam Co., Ltd.
958 Ikeuchi, Konan-cho, Takamatsu, Kagawa 761-4494 Japan

TUV Certified Management System
• ISO 9001, JIS Q 9001
• ISO 13485

REXXAM ELECTRONICS IRELAND CO. LIMITED
Dunmore Road, Drogheda County Louth, Ireland
Tel: +353-41-9839700 Fax: +353-41-9839702

SHIGIYA
SHIGIYA MACHINERY WORKS LTD.
GS Division <http://www.GrandSeiko.com>
5378 Minoshima, Fukuyama, Hiroshima, 721-8575 Japan
Tel: +81-84-981-5735 Fax: +81-84-953-6758

Printed in Japan 1005B2000N

REF/KERATO

AutoFocus AutoTrack Ref/Keratometer

AUTO REF/KERATOMETER
GR-3500KA
GR-3300K



Grand Seiko.com

3D Auto Measurement

With Auto Focus + Auto Track + Auto Start, the measurement is taken at the most proper position automatically without operating the joystick precisely nor pushing the Start Switch.
(GR-3500KA Only)



Auto Start

When the position and the focus are aligned, the measurement is taken automatically for the fixed times. Then, once both eyes are measured, print out will come out automatically.

Tilting Large LCD

Depending on the position of operator, the angle of LCD can be changed freely. With this Large Color LCD, the measurement can be done easily since the measurement information is shown in color and by icon.

Precise Measurement

With the newly designed optical system, the measurement accuracy has been improved.

User Friendly Fast Printer (with Auto Cutter)

For the replacement of paper, just put it in and close the cover. It is very simple and easy because there is no need to make any adjustment.



Simple Lock *1

Just turn the knob to stop the main unit temporary. It can be put on the sliding table safely.



Electric Chin Rest *2

Chin Rest moves electrically with the switch in front of operator.

Pupil Diameter Measurement

Simultaneously with refractive value, the pupil diameter at 0.1mm step (minimum 2mm) is shown on the LCD and can be printed out. It can measure the pupil diameter necessary for Multi Focus Contact Lens, IOL, and so on.

Advanced IOL Mode

After redesigning the whole system and making the improvement, even the cataract and the IOL implanted subject can be measured much more than before.

Data Output

With RS232C and USB Interface, the measurement data can be sent to the external devices. With the use of Data Collection Software (option), it is sent as spreadsheet data.

PD Measurement

Measures the Far PD and the Near PD is calculated accordingly (select working distance from 30, 35, 40, 45cm). The result can be printed out as NPD as well.


Small Pupil Measurement

Compare to the previous model, it measures smaller pupil of 2.2mm diameter to allow the measurement of elderly and the person with long eyebrows easier.

| Right eye data | | Date & time | |
|--|-------|-------------|---------|
| Message Area | | | |
| No. 00001 | | | |
| NAME 2009 05 28 13:10 | | | |
| VD=12 | | | |
| <R> | SPH | CYL | AX PUP. |
| | -4.75 | -0.25 | 90 5.0 |
| | -4.75 | -0.25 | 90 5.0 |
| | -4.75 | -0.25 | 90 5.0 |
| <R> | mm | D | AX |
| R1 | 7.59 | 44.50 | 120 |
| R2 | 7.57 | 44.50 | 30 |
| AVE | 7.58 | 44.50 | |
| | CYL | | |
| | 0.00 | | |
| RESI | AST | -0.25 | 90 |
| S | 7.84 | T | 7.66 |
| V | 7.51 | H | 7.55 |
| I | 7.84 | N | 7.71 |
| e (VER) | 0.554 | | |
| e (HOR) | 0.396 | | |
| e (AVE) | 0.475 | | |
| <L> | SPH | CYL | AX PUP. |
| | -4.50 | -0.75 | 90 5.0 |
| | -4.50 | -0.75 | 90 5.0 |
| | -4.50 | -0.75 | 90 5.0 |
| <L> | mm | D | AX |
| R1 | 7.59 | 44.50 | 120 |
| R2 | 7.57 | 44.50 | 30 |
| AVE | 7.58 | 44.50 | |
| Left eye data | | | |
| Average Eccentricity | | | |
| Horizontal Eccentricity | | | |
| Vertical Eccentricity | | | |
| Representative value (Displayed when the measurement is taken 3 times) | | | |

Corneal Peripheral Measurement

With the measurement of Corneal Peripheral, the base curve of contact lens can be chosen more precisely and the fitting of it become smoothly.



S(Superior) T(Temporal)
V(Vertical) H(Horizontal)
I(Inferior) N(Nasal)

Contact Lens Selection

With the addition of corneal peripheral, residual astigmatism, and pupil diameter data to the general refraction and corneal data, the accurate selection even for the progressive and the large size contact lens now become much easier.

Tonic Accommodation Measurement

The fatigueness during the near distance work and its recovery can be judged. At the same time, reliability of far point measurement can be confirmed. Traditionally, it must be measured in a special room by showing an "EMPTY FIELD".

| | | NAME 2009 05 28 13:10 | | Ref data | |
|----------------|-----|-----------------------|-------|----------|-------|
| Right eye data | <R> | SE | -5.43 | SE | -5.43 |
| | | TA | -5.48 | TA | -5.48 |
| | | TA-SE | -0.05 | TA-SE | -0.05 |
| Left eye data | <L> | SE | -5.43 | SE | -5.43 |
| | | TA | -5.48 | TA | -5.48 |
| | | TA-SE | -0.05 | TA-SE | -0.05 |

Grand Seiko.com
GR-3500KA

LCD Display (GR-3500KA)

VD Measurement Step Display of Measuring Place Map

Indication of Right/Left eye

Reticle mark

Indication of measurable pupil diameter

The result of radius measurement

No. of measurement

2nd principal meridian radius

1st principal meridian radius

2nd principal meridian axis angle

Interpupillary Distance

No. of measurement

The result of refractive measurement

Sphere S -4.75

Cylinder C -0.25

Axis angle A 90

Pupil Diameter PUPIL 4.0

•R1 7.85

R2 7.74

AX 90

PD=64 (R30+L34)